

Instrument Report

Instrument Details:

Name: Perkin Elmer Lambda 18 UV/Vis spectrometer

Model: Lambda 18

S/N: 75109

Installation date: 11 September 1996

Contact Information:

Chuanmin Hu (hu@marine.usf.edu)

College of Marine Science

University of South Florida


St. Petersburg, FL

USA

Document version: Version 1.0 (February 14, 2013)

I) Introduction

The Lambda 18 is a double beam UV/Vis spectrophotometer from Perkin Elmer, packing pre-aligned tungsten and deuterium lamps. It has a wavelength range of 185-900nm with a fixed bandwidth of 1 nm and wavelength accuracy of $\pm 0.15\text{nm}$. The photometric accuracy is $\pm 0.002\text{ AU}$ at 0.5 AU which translates to $\pm 0.046\text{ m}^{-1}$ when using a 10cm pathlength cell.



II) Calibration/Maintenance

Wavelength accuracy for this instrument has been checked intermittently using a Holmium Oxide filter on a near-annual basis.

III) Sample collection

Gelbstoff absorption spectra, $a_g(\lambda)$, was measured on filtered seawater samples collected at sea. Seawater was first filtered through GF/F filters (nominal pore size = $0.7\mu\text{m}$) to remove larger particles that tend to clog smaller pore size filters. These GF/F filtrates were then filtered through pre-rinsed $0.2\mu\text{m}$ nylon membrane filters. Samples were stored in clean 4oz amber borosilicate bottles and refrigerated when samples were processed within a week of collection or frozen if longer storage times were required.

IV) **Data processing**

Samples were allowed to warm up slowly to room temperature prior to being scanned. Duplicate scans of gelbstoff absorbance spectra, $A_g(\lambda)$, were measured from 200-800nm using 10cm quartz cells with Milli Q water as a reference. Raw individual scans were smoothed (7-point boxcar average) and averaged prior to being converted to $a_g(\lambda)$ as follows

$$a_g(\lambda) = \frac{2.3 A_g(\lambda)}{l}$$

where l is the cell pathlength. The average absorption from 690-710nm was then subtracted from the entire spectra to correct for residual scattering.

Reporting Notation:

abs_ag = absorbance coefficient (raw duplicate scans smoothed and averaged)

ag = absorption coefficient (null point corrected)